


PRINTING


PRINTING FIT EVAL - OPTION 2 (CURRENT SAMPLE)

Note: Current Sample Report will only print **ONE** sample round.


- 1. Navigate to the Fit Evaluation window.
- 2. Select desired sample round from the sample selection menu.




3. Click on Reports > Current Sample Report. PDF of **ONLY** that sample round will generate.



4. A second window will pop up automatically generating PDF report for printing. Select the Print Icon to enter preview mode.



5. Follow computer settings to either save as PDF file **OR** print PDF.



2/23/2021